

MASK FOR MOUTH-TO-MOUTH REANIMATION




Ambu

• **34216 LIFEKEY MASK**
Pocket mask for mouth-to-mouth reanimation. Its special one-direction valve allows oxygen to flow straight to the patient and minimizes risks of cross-contamination. Single-use, sterile, made in atoxic transparent material, easy to use and carry. Suitable for adults and children.
minimum order: pack 10 pcs.




GIMA

• **34218 CPR MASK - pocket resuscitator**
The CPR pocket mask is a non-invasive, latex free mask with an integral non-rebreathing valve and efficient filter to prevent direct mouth contact with the patient's face when providing mouth-to-mouth resuscitation. It can be used on adult, child and infant. Supplied in a compact plastic box with a corrugated tube for an easier use.



• **34211 VENTO AIRWAYS - Poly-bag**
• **34213 VENTO AIRWAYS - Key-holder box**
Disposable mouth to mouth airway with a one way polycarbonate valve and non-toxic pvc protective foil to avoid any direct contact between patient and rescuer. Available in 2 versions: in poly-bag or key-holder box. Instructions in GB, IT.

FORCEPS




GIMA code	STAINLESS STEEL FORCEPS
34001	Magill Forceps - 15 cm
34002	Magill Forceps - 20 cm
34003	Magill Forceps - 25 cm
34005	Collin Forceps - 20 cm
26745	Peaen Forceps - 14 cm
26720	Peaen Forceps - 16 cm
26710	Peaen Forceps - 20 cm

MOUTH OPENERS



GIMA code	PLASTIC DEVICES
34273	Teeth screw driver
34272	Helicoidal mouth opener
34274	Pinocchio
34274	Tongue holding forceps

STYLETS




Intubation stylets:
- coated with low friction material
- colour coded for easy identification
- malleable enough
- can be easily lengthened or shortened
- sold in box of 25
- individually packed
- non sterile
- made in nylon

GIMA code	Size	Colour	Ø mm	Length mm
34397	infant	blue	2.0	306
34398	child	grey	2.6	425
34399	adult	black	3.3	465

Single non sterile pouch - box of 25 pouches

EMERGENCY SCISSORS

GIMA code	EMERGENCY SCISSORS
34127	Emergency Scissors - 23.5 cm x 10.8 cm
34128	Universal Scissors - 15 cm
34126	Universal Scissors - 19 cm
34125	Lister Bandage Scissors - 11 cm
34129	Lister Bandage Scissors - 15 cm
34130	Lister Bandage Scissors - 18 cm



287



All purpose, first aid scissors, small enough to fit in the pocket of an emergency operator, easily cuts bandages, leather and plaster casts. Plastic coloured finger-contour handles for maximum comfort. Autoclavable at 136°C



Autoclavable Stainless steel AISI 420

Weight: 540 g

34127

Supplied with transport bag
Made in Italy