

AUTO NASAL IRRIGATOR

• 28087 AUTO NASAL IRRIGATOR

Easy-to-use device designed for gentle and effective nasal cleansing. Intelligent dual mode operation allows to switch between adult and child mode. Anti-backflow design prevents cross-infection, ensuring every flush is clean and hygienic. Pulse or continuous flow mode for personalized comfort and treatment efficiency. Stable Pulse Rinsing Technology controls the frequency of water flow pulses to penetrate deep into the nasal cavity.

The atomized particles are fine, dense and uniform. The flushing process is gentle and smooth without choking or irritation. Easy to disassemble and clean structure. It also has a self-standing base. Supplied with 1 adult and 1 child silicone tips, charging cable and manual in GB, IT.



28087


TECHNICAL SPECIFICATIONS

Spray rate:	continuous mode: ≥ 15.8 ml/min pulse mode: ≥ 6.3 ml/min
Air pump:	DC 6V, max inflation pressure ≥ 93 kPa
Max. capacity:	20 ml
Power supply:	3.7V Li-battery 800 mAh
Charging time:	about 2.5 hours
Size - weight:	64x52x167 mm - 256.9 g
Protection:	IP22

AUTO NASAL ASPIRATOR

• 28086 AUTO NASAL ASPIRATOR

Portable and lightweight electric nasal aspirator, designed for infants and newborns, offers gentle, fast and safe nasal cavity cleaning with minimal discomfort. 360° cleaning with no blind spots. Dual suction mode one for gentle cleaning and one for stronger suction - delivers effective, controlled and safe mucus removal. Low suction mode is for newborns, ensuring gentle and comfortable use. A one-way air intake system prevents backflow. All parts detach smoothly and the collection cup can be rinsed. The simple one-button operation provides full control: suction stops immediately when the button is released. Low-noise operation. Supplied with 4 different soft silicone heads: funnel-shaped tip, gourd-shaped tip, wide extended tip and wide tip, 2xAA batteries and manual GB, IT.



28086

TECHNICAL SPECIFICATIONS

Suction power:	≥ 200 mmHg (low suction) ≥ 350 mmHg (high suction)
Power supply:	2 x 1.5V AA batteries
Battery life:	1000 cycles
Size - weight:	155x56x44 mm - 190 g

NASAL ASPIRATOR

• 28085 NASAL ASPIRATOR

This nasal aspirator is designed to be used together with a compressor nebulizer to generate sufficient suction pressure for removing mucus from the nasal cavity. Intended for use in a general household setting. Suitable for adults and children starting from 1 month of age and older who suffer from upper airway conditions, such as nasal congestion, colds, allergies or respiratory conditions such as asthma, bronchitis or bronchiolitis. We recommend using it with AC or DC brushless compressor nebulizers, including GIMA models 28123 and 25058.



28085

TECHNICAL SPECIFICATIONS

Air flow:	4 l/min
Suction pressure:	0.14 ~ 0.2 bar
Min. input air pressure:	≥ 12 PSI

AROS SUCTION UNIT AND AROS BATTERY SUCTION UNIT



• 28050 AROS PORTABLE SUCTION UNIT - 1 l

• 28051 AROS BATTERY SUCTION UNIT - 1 l

This lightweight and easy to carry suction unit is ideal for home and emergency use. Modern and portable design. It's equipped with a rechargeable lithium battery that delivers up to 60 minutes of use on a full charge (28051 only).

The suction unit features a dual jar hook system, allowing quick changeovers without interrupting the procedure. It also supports disposable liner bags (optional) to simplify cleanup and minimize the risk of cross-contamination.

- upgraded vacuum gauge and adjustment knob for accurate suction control
- built-in overflow protection system safeguards the internal pump
- real time battery indicator to monitor power levels (28051 only)

STANDARD ACCESSORIES

Autoclavable jar 1,000 ml	Hose	Filters
Suction catheter	Fuse	Power cord
User manual GB, FR, IT, ES, PT, DE NL, PL, SE, HU, RO, GR		Suction tube 2m

TECHNICAL SPECIFICATIONS

	28050	28051
Suction power:	(28 $\pm$ 4) l/min	(27 $\pm$ 4) l/min
Vacuum regulating range:	20 kPa to maximum vacuum	
Max vacuum:	(85 $\pm$ 5) kPa	
Noise level:	≤ 65 dB(A)	
Power supply:	powered by an AC outlet	rechargeable lithium battery, 2600 mAh
Charging time:	-	4 hours
Operating voltage:	230 V, 50Hz	100-240 V, 50/60Hz
Input power:	120 VA	150 VA
Size:	480x165x285 mm	480x189x285 mm
Weight:	approx. 3.9 kg	4 kg
Protection:	IP21, Class II	

GIMA code	ACCESSORIES AND SPARES
28053	Autoclavable jar 1 l - spare
28054	Disposable suction liner 1 l - optional - box of 60
28055	Carrying bag - optional (only for 28051)
28056	Wall bracket for ambulance - optional (only for 28051)
28229	Spare filter
25468	Silicone tube 6x10 mm
28252	Adaptor for catheters