



## Technical data

| Pupil distance       |                                                                  |
|----------------------|------------------------------------------------------------------|
| Range                | 35mm to 80mm                                                     |
| Resolution           | 1mm                                                              |
| Accuracy             | $\pm 1$ mm                                                       |
| Gaze                 | 0° to 20°                                                        |
| Measuring distance   | 1m $\pm$ 5cm                                                     |
| Time per measurement | ~1s                                                              |
| Fixation target      | Light flash, attractive sound                                    |
| Date interface       | Wi-Fi, USB                                                       |
| Printer Interface    | USB                                                              |
| Battery              | Rechargeable lithium batteries, 6 hours of duration, Replaceable |
| Size                 | 180mm x 130mm x 110mm                                            |
| Display              | 5 inch touch screen                                              |
| Weight               | 0.8KG                                                            |
| Optional accessories | Camera tripods, printer, etc                                     |