



BIPOLAR CABLES FOR MAJOR BRANDS



Quality highly flexible silicone bipolar cables 3 m long (10").
Autoclavable at 134°C for minimum 20 times.
Recommended EtO sterilization for max extended life.

Suitable for most generators, with connector for European or US
forceps. Multilingual box and manual: GB, FR, IT, ES, PT, DE, SE PL,
RO,GR, HR, Arabic

<p>30642 BIPOLAR CABLE ADAPTOR</p>	<p>Connector to forceps</p>	<p>Connector to electro-surgical unit</p>	<p>Connector to forceps</p>	<p>Connector to electro-surgical unit</p> <p>Pin distance 29 mm</p>						
	GIMA Code	30632	30618	30640	30637	30643	30618+ 30642	30617	30644	30639
Connector to forceps	EU	US	EU	US	EU	US	EU	US	EU	EU
Compatibility brand/model	GIMA: MB 240, MB 380 VALLEYLAB, CONMED, CODMAN, OLYMPUS		GIMA: MB 120 F, MB 200 F, MB 200 D, MB 200T, MB 300 D, MB 400 D, MB 300T, MB 400T LED: SURTRON 120, 200, 200T, 300, 300T, 400, 400T VALLEYLAB, EMC, DOLLEY		GIMA: MB 122, MB 160, MB 200, MB 202		GIMA: MB 80 D, MB 120 D, MB 160 D LED: SURTRON 80, 120, 160		MARTIN, BERCHTOLD, AESCULAP, WOLF	ERBE, SELECT, DOWNS, SIEMENS



DISPOSABLE BIPOLAR FORCEPS - VALLEYLAB CABLE

- 30300 DISPOSABLE JEWELER FORCEPS - 11.5 cm - tip 0.5 mm
 - 30302 DISPOSABLE STRAIGHT FORCEPS - 20 cm - tip 1.0 mm
 - 30304 DISPOSABLE STRAIGHT FORCEPS - 15 cm - angled tip 1.0 mm
 - 30306 DISPOSABLE BAYONET FORCEPS - 20 cm - tip 1.0 mm
- Disposable bipolar forceps are supplied with 3 m cable.
Valleylab plug (same connector and compatibility as 30640).
Sterilization by Ethylene Oxide (EtO).





BIPOLAR SCISSORS

- 30310 BIPOLAR SCISSORS - 18 cm (7") - curved
Needs bipolar cable 30315
- 30315 ETHICON SILICONE BIPOLAR CABLE - 3 m (10')
For bipolar scissors. Connector 2.5 mm
Compatibility: same models as 30632, 30640, 30617*, 30643* (see above).
*Need adaptor 30642



BIPOLAR ARTERY SEALER

- 30320 BIPOLAR ARTERY SEALER - 26 cm (10.5")
Needs bipolar cable 30325
- 30325 SILICONE BIPOLAR ARTERY SEALER CABLE - 3 m (10')
Connector 4 mm, Banana plug.
Compatibility: same models as 30632, 30640, 30617*, 30643* (see above).
*Need adaptor 30642



FOOTSWITCHES AND BOARDS

GIMA code	FOOTSWITCHES FOR 50, 80, 106, MB 80D, 120D, 120F, 122, 132, 160, 160D, 200, 200D, 200F, 200 T, 202
30577	MB waterproof footswitch - 2m cable
30578	MB waterproof footswitch - 3m cable
	ANTIEXPLOSION FOOTSWITCH FOR MB 120D, 120F, 122, 132, 160, 160D, 200, 200D, 200F, 202
30602	Double antiexplosion, water proof MB footswitch
	ANTIEXPLOSION FOOTSWITCHES FOR MB 240, 300D, 300T, 380, 400, 400D, 400T
30603	Double antiexplosion, water proof footswitch - MB240-380
30604	Double antiexplosion, water proof footswitch - MB250-400-300D-400D
30608	Double antiexplosion, water proof footswitch - MB300-400 TOUCH
	MAIN BOARD
30557	Main board - specify model and serial number of diatermo with order
30558	Output board for MB 240, 250, 380, 400 - specify as above



It allows selection from pedal of CUT and COAG functions
Antiexplosion and waterproof.