

## ROLLATORS

Wheels  
Ø 10 cm (4")

CE  
Foldable



- **27721 ROLLATOR**  
Size 54x57 cm with adjustable height 52-87 cm.  
Weight: 8 kg - Load: 100 kg.



- **43132 ROLLATOR WITH SEAT**  
Blue painted steel adjustable height 79-93 cm frame with padded seat.  
Size: 52x49 cm. Weight 6.4 kg - Load: 100 kg.



Foldable to save space

## WALKING FRAME




- **27722 WALKING FRAME**  
Made of lightweight aluminium tubing, non-slip end and padded handgrips. Size 47x60 cm foldable with adjustable height 78-96 cm.  
Weight: 2.4 kg - Load: 100 kg.
- **27794 RUBBER TIP Ø 28 mm** - grey for 27722-3 - bag of 5 pcs.

## WALKING FRAMES WITH WHEELS



- **27723 WALKING FRAME - 2 WHEELS**  
Made of lightweight aluminium tubing, with padded hand grips, 2 non-slip ends and 2 with wheels.  
Size 51x64 cm, foldable with adjustable height 79-97 cm.  
Weight: 3 kg - Load: 100 kg.
- **27794 RUBBER TIP Ø 28 mm** - grey for 27722-3 - bag of 5 pcs.




- **43129 WALKING FRAME - 2 INTERNAL WHEELS**  
Made of lightweight aluminium tubing with padded hand grips, 2 non-slip ends and 2 with wheels. Foldable.  
Size 51x58 cm  
Adjustable height 80-90 cm.  
Weight: 2.5 kg - Load: 100 kg.



- **43130 WALKING FRAME - 4 WHEELS**  
Made of lightweight aluminium tubing, 2 front swivel wheels, 2 rear wheels with stop, with foam grips.  
Size 55x65 cm foldable with adjustable height 83-93 cm.  
Weight: 3.5 kg - Load: 100 kg.

## FOLDING SHOULDER WALKER

- **43136 FOLDING SHOULDER WALKER - adjustable height**  
Blue painted steel, foldable shoulder walker with chromed and padded, height adjustable shoulder support. Supplied with 4 castors, 2 of them with brake. Optional padded seat and brachial support. Made in Italy.
- **43137 PADDED SEAT**
- **43138 BRACHIAL SUPPORT**




Seat height: 50 cm

## ROTATING SEAT CUSHION

- **43040 ROTATING SEAT CUSHION**  
Cushion to assist to turn easily in a car, from the bed or on the chair.  
Size: Ø 40 cm



Slippery nylon filled with PU foam



Bottom in non-slip PVC fabric