

SERIES CHAGALL



CHAGALL 3 NEO - Art. 240000 - 3 sections

Quick height adjustment within 34 seconds.
Electrically adjustable backrest


Multifunctional professional treatment table with electrically adjustable height and backrest-tilt control.

For clinical examination, treatment, physiotherapy, massages

Version shown with optional footbar switch for peripheral control and castors

- 2 joint, 3 section table
- 2 electric motors with a thrust force of **10000 N** and **6000 N** (backrest)
- Manual control switch
- Safety control key
- Gas spring assisted leg section
- Trendelenburg and kyphosis position
- Electric height adjustment range 49 Δ 94 cm
- Fire-retardant padding and covering class 1IM, washable
- Safe working load 250 kg
- Dimensions: 196 x 65 x 49 Δ 94 cm
- Medical device

See detailed features page 48



OPTIONS

