


FOLDABLE BABY SCALE



• **27266 SOEHNLE 8320 FOLDABLE BABY SCALE**

Foldable baby scale for easy transport. Folding mechanism with magnetic lock. Easy to open: remove the base and pull the handles apart. Taring function: pads, blankets placed under the baby are tarable. Steady weight display through filtering and hold function. Large LCD display easy to read. Switchable to kg/lb.

Manual in: DE, GB, FR, ES, IT, GR, CZ, Arabic.

• **27271 CARRYING BAG**



TECHNICAL SPECIFICATIONS	
Capacity:	15 kg
Sensitivity:	10 g
Dimensions:	395x315x105 mm (closed) - 740x315x53 mm (open)
Weight:	1.9 kg
Alkaline battery 2x1.5 V AAA - not included	Warranty: 2 years

HOME CARE BABY SCALE



• **27306 SOEHNLE 8352 MULTINA DIGITAL BABY AND CHILD SCALE**

Well-designed, easy to handle plastic baby scale. - suitable for children up to 20 kg removing baby tray - easy to read LCD-display with 13 mm figures - simple operation with one button - tare function to calibrate for baby blankets/textiles - fully-electronic scale: no wear and tear parts, no maintenance - hold function: weight value will be held for 10 seconds after unloading - automatic power-off

- **27307 TRAY FOR 8352 - spare**
- **27264 BABY HEIGHT ROD - optional**
- see page 87



TECHNICAL SPECIFICATIONS			
Capacity:	20 kg	Sensitivity:	10 g
Dimensions:	601x346x92 mm	Weight:	3.5 kg
Battery:	9 V	Warranty:	2 years

BABY AND CHILD SCALE



• **27312 BABY AND CHILD DIGITAL SCALE**

Multifunctional scale with removable light plastic tray. It can be used with tray for babies or without for children. Weight difference display function with memory: you can see the gain of weight of the baby after feeding. Precise scale with easy-to-read LCD display. Auto zero, auto tare, auto hold function. Switchable to kg/lb/st.



TECHNICAL SPECIFICATIONS			
Tray dimensions:	525x305x63-95 mm	Capacity:	20 kg
Pack dimension:	540x300x110 mm	Sensitivity:	5 g
Battery operated:	2xCR2032 lithium batteries	Weight:	2 kg

GIMA ELECTRONIC BABY SCALE



• **27313 GIMA ELECTRONIC BABY SCALE**

Handy, reliable and well-designed electronic baby scale. Provided with large LCD display 66x28 mm and kg/lb switch. Supplied with height tape meter and washable sheet. Multilingual user manual and gift box: GB, FR, IT, ES, PT, DE, Arabic



TECHNICAL SPECIFICATIONS			
Capacity:	20 kg / 44 lbs	Sensitivity:	10 g
Dimensions:	595x340x90 mm	Weight:	4 kg
Battery operated:	4x1.5 V AA	Auto-off function	

FAMILY BABY SCALE

• **27302 FAMILY BABY SCALE**

Baby scale for domiciliary and family use. Capacity 20 kg/44 lbs. Scale division: 50 g (0-10 kg), 100 g (10-20 kg) Plastic tray. Dialplate in kg and lbs. Manual zero adjustment. Multilingual manual: GB, FR, IT, ES, PT, DE, Arabic.

