

Hot/cold pack MD CE

The Moves Hot/Cold Packs are used for hot as well as cold therapy. This dual applicability makes them a strong alternative to ordinary heating pillows and ice packs. The versatile, reusable packs deliver therapy effectively and economically. The use of the Moves Hot/Cold Pack is unlimited and its material is non-toxic. Even below -20°C it remains soft. The Moves Hot/Cold Pack can be bent in any angle, allowing maximum coverage of joint and muscle. All Hot/Cold Packs can be heated (microwave) and can be frozen (freezer). The temperature holds for approximately 30 minutes.

CLASSIC - GEL

The classic (and economical) transparent Hot/Cold Pack. The special gel is non-toxic and biodegradable.

✓ Sold per 10 pcs

🏠 For home use only!



STANDARD - GEL

Our standard Hot/Cold Pack is a nylon covered gelpack. The special gel is non-toxic and biodegradable.

✓ Sold per 10 pcs



SOFT TOUCH - CLAY

The soft touch provides a user-friendly and frost-free surface that allows usage in **direct contact with the skin**. Filled with clay, a 100% natural product.


✓ Sold per 10 pcs



PROTECTION SLEEVE

100% polyester Protection Sleeve for Hot/Cold Packs. Machine washable. Available in 3 sizes (S | M | L).

✓ Sold per 10 pcs



TRI-SECTIONAL HOT/COLD PACK

The soft touch tri-sectional Hot/Cold Pack treats multiple zones at once (neck and shoulders).

✓ Sold individually & per 10 pcs



XXL HOT/COLD PACK

Ideal for pain relief on shoulders and back. This Hot/Cold Pack has a nylon cover.

✓ Sold individually



Dimension		Moves Hot/Cold Packs				Unit
		Classic (Gel)	Standard (Gel)	Soft Touch (Clay)	Protection Sleeve	
Small	15 x 25 cm	07-010102	07-010202	07-010302	07-020102	10 pcs
Medium	20 x 30 cm	07-010103	07-010203	07-010303	07-020103	10 pcs
Large	25 x 35 cm	07-010104	07-010204	07-010304	07-020104	10 pcs
XXL	33 x 47 cm	n/a	07-010206	n/a	n/a	1 pc
Tri-Sectional	40 x 20 cm	n/a	07-010211	07-010311	n/a	10 pcs
Tri-Sectional	40 x 20 cm	n/a	n/a	07-010312	n/a	1 pc