

20D




PRIMARY APPLICATION Industry Standard General Diagnostic Lens

- + Perfect balance of magnification and field of view makes this lens the most popular choice for general diagnostic exams.
- + Dynamic examination allows viewing of the peripheral retina while a primary position gaze enables a central retinal exam
- + Also available in autoclave sterilizable (ACS*) design (see page 53) or single-use design (see page 58)

46°/60°
FIELD OF
VIEW

3.13x
IMAGE
MAG

0.32x
LASER
SPOT MAG



Lentila diagnostic 20D Volk® echilibru perfect între mărime și precizia câmpului vizual în evaluarea și diagnosticul oftalmologic general. Lentile de diagnostic Volk® 20D BIO optica de calitate superioară, dubla asferică, considerată cea mai recunoscută lentilă în salile de diagnostic oftalmologic modern.

Lentila diagnostic 20D BIO oferă imagini detaliate, câmp vizual excelent, experiență confortabilă pentru medicul oftalmolog în manipularea lentilei, este ideală pentru examinările centrale ale retinei de nivel înalt. Lentilele Volk® reprezintă alegerea ideală pentru profesioniștii din domeniul oftalmologiei în evaluarea și diagnosticul afecțiunilor oculare.

Lentile diagnostic: Volk® 20D BIO

Câmp vizual: 46° / 60°

Mărime imagini: 3.13x

Mărime spot laser: 0.32x

Distanța de lucru: 50 mm

Culori disponibile: negru, la cerere - albastru, auriu, verde, violet, roșu, argintiu

