

60D



V60C

PRIMARY APPLICATION


High Magnification View of the Posterior Pole

- + High magnification lens for detailed optic disc and macula imaging
- + High magnification enables detection of small defects and subtle changes in retinal abnormalities
- + Ideal diameter for use in the orbital area
- + Dilation is required to obtain optimum retinal imaging

68°/81°
FIELD OF
VIEW

1.15x
IMAGE
MAG

0.87x
LASER
SPOT MAG



Lentila diagnostic 60D, ofera imagini clare, detaliate, destinata pentru evaluarea si diagnosticul oftalmologic in tratamentul afectiunilor oculare, anomaliiilor retiniene, degenerescenta maculara legata de varsta (DMLV, este o afectiune oculara care se manifesta prin pierderea acuitatii vederii centrale), urmand raportul pahar-disc la pacienti si hemoragiile capilare.

Lentile de diagnostic Volk® optica de calitate superioara, design dublu asferic, imagini centrale detaliate ale retinei, dilatarea pupilei este necesara pentru o imagine optima. Lentila diagnostic 60D Volk ofera marire mare pentru imagini detaliate ale discului optic si ale maculei pentru a detecta mici defecte si modificari subtile ale anomaliiilor retinei. Esentiala pentru oftalmologi în diagnosticul oftalmologic si tratamentul afectiunilor oculare, diametru perfect pentru utilizare în zona orbitala.

Lentile diagnostic: Volk® 60D

Camp vizual: 68° / 81°

Marire imagini: 1,15x

Marire spot laser: 0,87x

Culori disponibile: negru

