

78D



PRIMARY APPLICATION


General Diagnosis and Laser Treatment

- + Ideal balance of magnification and field of view
- + Optimally designed for use within range of motion of all slit lamps
- + Offers clear and large views of the central mid-retinal regions
- + Dilation is required to obtain optimum retinal imaging
- + Ideal general lens for doctors who regularly cater to populations prone to glaucoma and other posterior pole abnormalities

81°/97°
FIELD OF VIEW

0.93x
IMAGE MAG

1.08x
LASER SPOT MAG



Lentile de diagnostic 78D Volk®, ofera imagini detaliate si echilibru excelent de marire a campului vizual, pentru diagnostic oftalmologic modern. Lentila diagnostic 90D Volk®, optica de calitate superioara, design dublu asferic, ofera imagini detaliate mari, clare ale regiunilor centrale medii retiniene si alte anomalii ale polului posterior.

Lentila 78D Volk® recomandata medicilor oftalmologi pentru diagnostic oftalmologic general, esentiala în evaluarea si diagnosticul afectiunilor oculare, glaucom etc. Lentilele Volk® reprezinta alegerea ideala pentru profesionistii din domeniul oftalmologiei.

Lentile diagnostic: Volk® 78D

Camp vizual: 81° / 97°

Marire imagini: 0.93x

Marire spot laser: 1.08x

Distanța de lucru: 8 mm

Culori disponibile: negru, albastru, auriu, verde, violet, rosu, argintiu

