

3-Mirror



No Flange:
V3MIR (shown)

ANF+ Flange:
V3MIRANF+

PRIMARY APPLICATION


3-Mirror, Acrylic Design for Anterior, Peripheral, and Equatorial Viewing (Goldmann-style Lens)

- + 3-mirror design provides the same anterior chamber angle, central, equatorial, and peripheral retinal views as our G-3 Gonio lenses, but in a light-weight acrylic design while still providing Volk quality optics
- + Advanced no fluid (ANF+) flange only requires a coupling fluid during laser procedures
- + Not recommended for SLT. We recommend a Volk Rapid SLT® or SLT lens instead (page 39)

60°/66°/76°
MIRROR
ANGLES

0.90x
IMAGE
MAG

1.11x
LASER
SPOT MAG



Lentila diagnostic Volk® 3 Mirror, acrilica, usora, design compact si claritate optica superioara, lentila Volk® trei oglinzi, asigura o vizualizarea clara si detaliata a unghiului iridocornean, a retinei central periferice si a retinei medii periferice/periferice.

Lentile de diagnostic Volk® 3 Mirror ofera imagini detaliate si unghiuri multiple pentru evaluarea ochiului. Lentila Volk® asigura diagnosticul precis al patologieilor oculare, esentiala pentru oftalmologi în examinarile detaliate. Lentila diagnostic Volk® 3 Mirror este un instrument de încredere în diagnosticul oftalmologic si tratamentul afectiunilor oculare.

Lentile diagnostic: Volk® 3 Mirror

Camp vizual: 60°/66°/76°

Marire imagini: 0.90x

Marire spot laser: 1.11x

Distanța de lucru: 3 mm / 18 mm (ANF+)

Culori disponibile: albastru

