

QuadrAspheric®



Flange:
VQFL (shown)

No Flange:
VQFLNF

ANF+ Flange:
VQFLANF+

PRIMARY APPLICATION | SMALL PUPIL PRP

Wide Field View for Pan Retinal Examination and Laser Treatments in Small Pupils

- + High resolution imaging of the peripheral retina with small pupil capability for patients who do not accommodate dilation well
- + Excellent general diagnostic and laser treatment lens
- + Flanged version recommended for laser for optimal stability on cornea

120°/144°


FIELD OF
VIEW

0.51x

IMAGE
MAG

1.97x

LASER
SPOT MAG



Lentila diagnostic QuadrAspheric® solutia suprema, pentru profesionistii din domeniul oftalmologiei în evaluarea si diagnosticul oftalmologic al afectiunilor oculare, retiniene periferice.

Lentile de diagnostic Volk® designul sau unic permite vizualizarea campului larg pentru evaluarea si tratarea pacientilor cu risc ridicat, predispusi la închiderea unghiului si neovascularizarea irisului, pupile mici si unde dilatarea trebuie evitata. Lentila diagnostic QuadrAspheric® Volk optica de calitate superioara, echipat cu flansa mare, ofera stabilitate si control de neegalat în timpul procedurilor cu laser.

Lentile diagnostic: Volk® QuadrAspheric

Camp vizual: 120° / 144°

Marire imagini: 0.51x

Marire spot laser: 1.97x

Culori disponibile: negru

