

PICASSO QUATTRO

Electric osteopathy table with four sections. The height is variable and assisted by an electric motor. Recommended also for massages and medical exams.



Medical
Device



Dimensions
198 x 64 x 50 h92 cm



Adjustable height
da cm 50 cm min a 92 cm max



Safe working load
140 Kg



Maximum capacity
180 Kg



Unimpeded passage
for patient lifting



The photo depicts the product complete with optional accessories

FEATURES

Gas piston assisted adjustment of the backrest and footrest

Electric motor with a thrust force of 8000 N

Push-button control panel

Safety control key

Class IM1 fire-retardant padding and upholstery, washable
