

PICASSO COMFORT

Two-section physiotherapy electric table with forward and backward tilting of the backrest, and electrically assisted height adjustment. Recommended for therapy sessions, massages and medical exams. Equipped with paper roll holder.



Medical
Device



Dimensions
209 x 64 x 50 h92 cm



Adjustable height
da cm 50 cm min a 92 cm max



Safe working load
160 Kg



Maximum capacity
200 Kg



Unimpeded passage
for patient lifting



Version presented with perimeter command and wheel options

FEATURES

Two-section table with with 84 cm forward tilting of the backrest only, assisted by gas piston.

Electric motor with a thrust force of 8000 N

Push-button control panel

Safety control key

Class 1IM fire-retardant padding and upholstery
