



MAIN
MEDITECH

SVF-7000
Main Perimeter

Model Parameter	SVF-6000	SVF-7000
DB Value Range displayed	0~40db	
Stimulus Intensity Range	1asb~10000asb	
Stimulus Size	Goldmann III	
Stimulus Color	White	
Stimulus Duration Time	200ms/ Adaptive	
Max Measurement range	90°	
Test Distance	30cm	
Test Mode	Static	
Stimulus quantity	456 points	
Stimulus Plane	Aspherical surface	
Background Illumination	White 31.5 asb (10cd/m2)	
Customize test program	/	●
Auto eye recognition	/	●
Blindzone range test	/	●
Stimulus Source	LED	
Chinrest Control	Touch Screen,Keyboard and mouse control,Button control	
Head Movement	Chinrest and forehead-rest move simultaneously	
PC Configuration	CPU:1.6GHz; RAM:4GB; Hard disk : 120GB; Operation system: Windows 10; LCD: 14 " Capacitive Touch screen; Networking: Ethernet&WIFI;	
Dimension	480mm(L)*430mm(W)*528mm(H), 14kgs	
Packaging	550mm(L)*480mm(W)*640mm(H), 19kgs	
Others	Input Voltage: AC 100 - 240V,50 ~ 60Hz	