

STAINLESS STEEL HOLLOWARE - AUTOCLAVABLE 121°C - AISI 304



STAINLESS STEEL HOLLOWARE
Stainless steel AISI 304 holloware with minimum thickness 0.6 mm. Autoclavable at 121°C.

GIMA code	INSTRUMENT TRAYS - ANGLED EDGES
	Without cover
26603	306x196x h 50 mm
26604	355x254x h 50 mm
26597	380x304x h 50 mm
26598	440x320x h 64 mm
GIMA code	With cover
	26607 306x196x h 50 mm
	26608 355x254x h 50 mm
	26609 440x320x h 64 mm



GIMA code	INSTRUMENT TRAYS ROUND EDGES
	Without cover
26599	223x126x h 45 mm
26602	264x172x h 47 mm
GIMA code	With cover
	26605 223x126x h 45 mm
	26606 264x172x h 47 mm

GIMA code	DENTAL TRAY
26600	208x109x h 15 mm
INSTRUMENT TRAYS	
26601	210x160x h 25 mm
26629	313x218x h 31 mm

GIMA code	MAYO TRAYS - SHALLOW
26613	254x165x h 18 mm
26614	350x252x h 16 mm
26615	432x295x h 19 mm



- 26592 BED PAN ROUND WITH COVER - 320x85 mm



- DRESSING JAR WITH LID
- 26623 DRESSING JAR - 0.5 l - Ø 102x65 mm
- 26624 DRESSING JAR - 1 l - Ø 103x126 mm
- 26625 DRESSING JAR - 2 l - Ø 127x161 mm



- 26626 FORCEPS JAR - Ø 55xh 140 mm
- 26627 FORCEPS JAR - Ø 55xh 180 mm
- 26628 THERMOMETER JAR - Ø 33xh 80 mm



- 26590 BED PAN FRACTURE TYPE - with ring - 418x292x85 mm



- 26594 BED PAN FRACTURE TYPE - with tube handle - 418x292x85 mm



- 26593 BED PAN PERFECTION TYPE 354x289x97 mm



- 26578 MALE URINAL - 1.2 l



- 26579 FEMALE URINAL - 1.2 l



- 26591 MALE VERTICAL URINAL Ø 114x220 mm - 1.5 l