



Product Code	34350
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B
EAN/UPC	6945040100232

Description

GIMA WRIST PULSE OXIMETER - with software

Small and lightweight wrist pulse oximeter for continuous monitoring of SpO2 and pulse rate, suitable for sports monitoring or at home, hospital and healthcare facilities.

- display of SpO2 value - display of pulse rate and bar graph
- display of pulse waveform - low battery indication
- adjustable screen brightness - PR sound indication
- data memory function - clock function and charging function
- data can be uploaded to the terminal equipment by wired mode
- sound prompt function for over limit, probe-off or finger-out and low battery

Supplied with adult SpO2 probe, USB cable, software multilingual manual in GB, FR, IT, ES, DE, NL, SE, PL.

Technical Specifications

SpO2 range: 0~100% (resolution: 1%) / accuracy: ± 2% when range is 70~100%; below 70% unspecified
Pulse rate range: 30~250 bpm (resolution: 1 bpm) / accuracy: ± 2 bpm or ± 2 %
Wavelength: red light: 660 nm; infrared light: 905 nm
Memory: 99 group of data under auto mode, 72 hours total up to 24 hour data under manual mode
Battery: 3.7 V rechargeable lithium battery
Battery life: > 500 times
Size - weight: 61x56xh 24 mm - 60 g
Protection: type BF applied part, IP22