

FOLDABLE COMMODE WHEELCHAIR



Product Code	43202
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2693045
CND	Y091212
EAN/UPC	8023279432022
PARAF	950833703

Description

FOLDING COMMODE WHEELCHAIR

- one-button foldable commode
- EVA seat and backrest with soft surface and odour free
- detachable seat and backrest pads
- rust-proof frame, made of aluminum alloy
- equipped with removable toilet pan

Packing size: 87x23.5x63 cm

Max. load: 125 kg

Weight: 9.4 kg

Multilingual manual: GB, FR, IT, ES, PT, DE, SE, PL, HR, RO, GR.

PU rear wheels, Ø 30.5 cm

Front castors with brake, Ø 10 cm