

# IV Pole Normal care

Material: chromed steel

I-N41221

- + effective infusion care by hanging safely IV-bottles and bags
- + mobile use: at the right time at the right place



## Product description

Material:	chromed steel
Colour of plastic parts:	blue
Bottle holder:	plastic 4 safety bottle hooks - 2 kg / hook
Adjustment:	screw safety height adjustment
Regulating range:	1.350 - 2.150 mm
Tube diameter:	25 / 18 mm
Length of outer tube:	1.000 mm
Base diameter:	635 mm
Base leg profile:	25 x 25 mm
Castors:	twin castors Ø 50 mm, according to DIN EN 60601 electrically conductive and lockable
Base construction:	weighted
Total- / base weight (net):	appr. 6,1 kg / 4,5 kg



# IV Pole Normal care

Material: chromed steel



## Application

- + optimal hygienic characteristics (wipe disinfection)
- + very easy and fast assembly through a plug in and out system with only 1 screw connection
- + comfortable height adjustment through a screw-height adjustment with an integrated, surface-gently mold clamping element

## Technology

- + all plastic parts are light-resisting
- + castors with non-marking tread
- + integrated pull-out guard avoids separation of the adjustable tubes

## Economy

- + complete recyclability

## Safety

- + enough safe distance in between both tubes once the pole reached its lowest height
- + an integrated base-weights facilitates a safer attachment of various pumps and other medical equipment
- + no risk of injury from the rounded shape of the hook