


Height rod foldable 5003



- ▶ Handy height rod for mobile use
- ▶ For measuring yourself and others

Characteristics

- ▶ Ultrasound measuring technology, temperature-compensated
- ▶ Tilt sensor for correction if measuring instrument is not horizontal, resulting in an error message
- ▶ Low space requirement due to fold function
- ▶ Plastic housing
- ▶ LCD display, 9 mm digit height
- ▶ Acoustic acknowledgement at end of measurement and reading displayed
- ▶ Units in cm or inch
- ▶ Battery operation with 2 x 1,5 V AA batteries (included in scope of delivery)

Order No.	Dimensions mm	Measuring range cm	Division cm	Own weight g
5003				
5003.01.001	345 x 92 x 28	50 – 240	1	250