

GIMA A202 EAR INFRARED THERMOMETER
GB,ES,PT,NL,SE,LT,PL,CZ,HU,HR,RO,GR



Product Code	25592
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

Description

GIMA A202 EAR INFRARED THERMOMETER GB, ES, PT, NL, SE, LT, PL, CZ, HU, HR, RO, GR

Portable and easy-to-use infrared thermometer supporting both ear and environment modes. It features a high-quality sensor providing accuracy within 0.2°C, delivering results in just one second. To ensure reliable monitoring and tailored results, an advanced algorithm provides specific measurement modes for infants aged 0 to 36 months and users aged 36+ months.

- three-colour backlit LCD display (green, yellow, red) for immediate visual feedback
- 25-reading memory for easy tracking
- switch between °C and °F
- audible reminder
- low battery indicator

Supplied with 20 probe covers and 2xAAA batteries and multilingual manual.

Technical Specifications

Measuring range:
body: 32-43°C (89.6-109.4°F)
environment: 0-100°C (32-212°F)

Accuracy:
body: ± 0.2°C (0.4°F)
environment: ± 1°C (33.8°F)

Sound reminder: 38-43°C (100.4-109.4°F)
Size - weight: 143x38x49 mm - 78 g