



The familiar way



**Compactly combined instruments
that pay dividends.**

Doctor's bags and cases

med-kit® I

The med-kit® I doctor's case which contains high quality Riester products is an indispensable piece of basic equipment for every physician. There is a choice of the palm-style blood pressure monitoring device ri-san® which is suitable for right or left-handed individuals or the precisa® N, an ideal mobile device. The set also includes the very light duplex® aluminum stethoscope with particularly good acoustic performance, the high-quality buck percussion hammer for exact testing of reflexes and the sturdy and convenient sized diagnostic penlight fortelux® N.



Order numbers med-kit® I:

- Aneroid sphygmomanometer precisa® N plastic/metal
- Stethoscope duplex® aluminum
- Diagnostic penlight fortelux® N
- buck percussion hammer
- Dimensions:
L 31.5 cm x W 29 cm x H 10.5 cm
- No. 5280 black (see illustration)**

- Aneroid sphygmomanometer ri-san®
- Stethoscope duplex® aluminum
- Diagnostic penlight fortelux® N
- buck percussion hammer
- Dimensions:
L 31.5 cm x W 29 cm x H 10.5 cm
- No. 5280-02 slate grey**
- No. 5280-03 blue**
- No. 5280-04 saffron**
- No. 5280-05 green**

Replacement parts included:

- precisa® N/ri-san® protective glass
- duplex® ear tips
- duplex® diaphragm
- fortelux® N light
- buck rubber tips



duplex®



Tube-screw connection for quick exchange of cuffs with the precisa® N.



Convenient regulation of air-release valve by simple finger-tip control with the ri-san®.



precisa® N



fortelux® N




ri-san®



buck

med-kit® II

With the quality doctor's case med-kit® II you have numerous reliable and durable instruments to hand, which are needed for mobile use. Alternative equipment includes, amongst others, the uni® I otoscope 2.7 V in a sturdy metal design with vacuum illumination and adjustable light intensity using a rheostat or the ri-mini® F.O. otoscope 2.5 V with halogen illumination and fiber optics.



Order numbers med-kit® II:

- Aneroid sphygmomanometer precisa® N plastic/metal
- Otoscope uni® I
- Stethoscope duplex® aluminum
- Diagnostic penlight fortelux® N
- buck percussion hammer
- Dimensions:
L 38 cm x W 30 cm x H 8 cm
- No. 5281 black (see illustration)**

- Aneroid sphygmomanometer ri-san®
- ri-mini® F.O. otoscope
- Stethoscope duplex® aluminum
- Diagnostic penlight fortelux® N
- buck percussion hammer
- Dimensions:
L 38 cm x B 30 cm x H 8 cm
- No. 5281-02 slate grey**
- No. 5281-03 blue**
- No. 5281-04 saffron**
- No. 5281-05 green**

Replacement parts included:

- precisa® N/ri-san® protective glass
- uni® I/ri-mini® ear specula in various sizes
- uni® I light
- duplex® ear tips
- duplex® diaphragm
- fortelux® N light
- buck rubber tips



buck



precisa® N



ri-san®



fortelux® N



uni® I



ri-mini®



duplex®

The first-class equipped doctor's case med-kit® III offers a full range of high-quality Riester quality products. With this useful companion all important instruments will always be to hand. In the med-kit® III you have, amongst others, the high quality F.O. otoscope and ophthalmoscope instrument heads from the excellent Riester ri-scope® line with particularly bright 2.5 V HL halogen illumination and the robust battery handle type C with rheostat for regulating light intensity. The high performance stethoscope duplex® de luxe is convincing with its excellent acoustics for precise auscultation of adults and children.




Order number med-kit® III:

- Aneroid sphygmomanometer precisa® N plastic/metal
 - ri-scope® F.O. otoscope
 - ri-scope® ophthalmoscope
 - Battery handle type C
 - Stethoscope duplex® de luxe® stainless steel
 - Diagnostic penlight fortelux® N
 - buck percussion hammer
 - Aluminum tuning fork c 128 Hz acc. to Hartmann with clamps
 - Dimensions: L 38 cm x W 30 cm x H 8 cm
- No. 5282**

Replacement parts included:

- precisa® N protective glass
- ri-scope® F.O. otoscope fast plug-in connector
- ri-scope® lights
- ri-scope® F.O. otoscope ear specula in various sizes
- duplex® de luxe ear tips
- duplex® de luxe diaphragms
- fortelux® N light
- buck rubber tips



D+	1	2	3	4	6	10	15	20
D-	1	2	3	4	6	10	15	20



ri-scope® F.O. otoscope



ri-scope® ophthalmoscope



Battery handle type C



duplex® de luxe



Tuning fork c 128 Hz



buck

A detailed description of the med-kit® products can be found in the corresponding chapter of our product catalogue.

Rudolf Riester GmbH & Co. KG
P.O. Box 35 • DE-72417 Jungingen
Germany
Tel.: (00 49) 74 77/92 70-0
Fax: (00 49) 74 77/92 70-70
info@riester.de
www.riester.de